

Sound LOGIQ™ e Vet Ultrasound

Simple, Faster and even more Precise.



Sound LOGIQ™ e

Portable imaging taken to the next level

The new LOGIQ e brings a wave of innovation to our most versatile and robust laptop ultrasound system.

Advanced Hardware

- New intuitive user interface with fewer controls
- Integrated carrying handle
- Larger high-resolution viewing monitor
- Electronic TGC's
- 1.6 lbs. lighter than the previous generation product
- Up to 90 minutes of battery life
- Whisper quiet cooling fans

Accelerated Software

- Plug and Play veterinary presets for out of the box, high-end imaging
- Needle Guidance software
- Patient Follow Up Tool with Fusion for serial studies
- Auto EF for quick Ejection Fraction measurements
- Strain Elastography

New Probe Technology

A new L4-20t-RS XDclear transducer has been added to the LOGIQ e. This transducer represents a leap forward in image quality and usability with remote operation via configurable buttons. Added flexibility and functionality for an ever-changing clinical environment.



Features

- Supports XDclear transducers
- Animal patient ID screens
- Optimized veterinary presets
- Continuous Wave Doppler
- Customizable Scan Assistant protocol software
- Patient Follow Up Tool for serial studies
- LOGIQ VIEW
- Needle Guidance software
- Probe check software




ML6-15-RS

Frequency range from 5-15MHz




12S-RS

Frequency range from 4.1-12MHz




Elastography

Provides an estimate of strain or tissue deformation to identify pathology sooner.



B-Steer/Needle Guidance


Enhances the ability to see both the target and needle with independent needle and B-mode gain controls.



Follow Up/Fusion

Optimize treatment with serial studies. Align anatomy from previous studies during live imaging while synchronizing all imaging parameters.

GE Exclusive



Hi-Res PDI

Diagnose and follow treatment of inflammation by detecting and quantifying the presence of microvasculature and blood flow.

GE Exclusive

Simple

Your work can be complicated so we made LOGIQ e simpler. LOGIQ e's specialized software and transducers help you keep up with treatment advances and simplify interventions. The updated design is both ergonomic and functional for easy connectivity and operation.

Faster

We've made scanning Faster, giving you more time for your patients. Get up to speed quickly with everything you need directly on the keyboard. Scan Assistant makes quick work of tedious tasks to keep you focused on providing the best care. Auto EF for accurate and instant ejection fraction results.

Even More Precise

Image quality matters. So we made LOGIQ e even more precise. As image quality continually improves, so does the precision of diagnostic, guidance, and monitoring treatments. That's why LOGIQ e incorporates advanced imaging technology to help reveal the whole story.

8C-RS



12L-RS



L4-20t-RS



6S-RS



10-22-RS



4C-RS



C1-5-RS



3Sc-RS



L8-18i-RS



LK760-RS



